




1 - Some tools of the trade

Lab 1.1: Using a micropipette



Check reagent (RD) → Review use of micropipette → Practise setting micropipette volumes at 2.0 μL , 5.5 μL , 12.4 μL and 20.0 μL



Micropipette Practice:					
Name	20.0 μL	15.0 μL	10.0 μL	5.0 μL	2.0 μL
1. _____	☐	☐	☐	☐	☐
2. _____	☐	☐	☐	☐	☐
3. _____	☐	☐	☐	☐	☐
4. _____	☐	☐	☐	☐	☐

Micropipette the volumes onto your laminated sheet as accurately as possible